



Film and Paper Digitizer



Specifications

Resolution	Up to 600 dpi
Sensor	High Pixel CCD
Optical Density (DMAX)	4.7
Bit Depth	16 bits, 65536 shades of gray output
Film / Paper Size	Min: 1.96" x 2.17" (50 mm x 55 mm) Max: 14" x 35" (355.6 mm x 889 mm)
Scan Rate	11.5 ms/line
Interface	Hi-Speed USB 2.0
Software	Windows XP and Vista interface compatible with Omni-Vue™
Voltage Requirements	AC 100V to 240V, 47-63 Hz, 1.5A Max (Input)
Power Consumption	54.9 Watts (Max)
Operating Environment	50° – 104° F (10° – 40° C)
Relative Humidity	20% to 85%
Light Source	Cold Cathode lamp
Lamp Life	>= 50,000 hours
MTBF	>= 50,000 hours
Weight	19.82 lbs. (9 kg)
Dimensions (L x W x H)	10.2" x 18.7" x 9.3" (259 mm x 475 mm x 236 mm)
Certifications	FDA, CE, FCC Class B, ISO 9001, ISO 14001

Specifications are subject to change without notice