



Managed Services

UltraRAD Corporation offers an innovative solution to outsource all of your digital imaging needs. Our team has over 20 years of experience and hands-on implementation and integration of Teleradiology, PACS, and medical informatics solutions to provide your facility with both comprehensive solutions and support.

Offsite Archiving

UltraRAD provides a scalable, cost-efficient storage solution that will enable fast and accurate access to critical patient medical data. Our team is experienced in assisting organizations in their move towards a filmless and protected environment.

Disaster Recovery

UltraRAD provides a scalable, cost-effective storage solution that enables access to patient data when you need it most. With our Disaster Recovery Service our customers have the peace of mind of knowing their critical data is maintained and secure for their access when disaster strikes.



Hosted Web Services

UltraRAD Hosted Web Service is the perfect solution for allowing unlimited access to your critical patient information. UltraRAD provides a scalable, cost-effective web based viewing solution that enables access to patient data when you need it most.

Our Solution delivers medical image and information viewing and distribution without the high capital cost, maintenance fees, technical upgrades and staffing support typically associated with traditional PACS systems. You can easily expand the system to include our Off-Site Archiving Service or Disaster Recovery Service.

The UltraRAD platform allows for easy integration with any DICOM-compliant system. Our solution scales from a single modality to multi-site PACS installations.

Integrated Solutions for Medical Information Management



Integrated Solutions for Medical Information Management

UltraRAD Corporation has been serving the medical community for over 20 years. We have built a dependable reputation in the medical information and image management industry by developing and delivering innovative technology solutions that facilitate the medical community's ability to manage their business. As the demand for information management within healthcare grows, UltraRAD understands that workflow and interoperability between disparate sources is vital. UltraRAD has developed applications that allow each potential customer from any size facility the best integrated workflow in digital imaging. UltraRAD's full suite of solutions can be utilized to serve your specific workflow requirements.

Information Acquisition and Management



UltraBROKER

UltraBROKER is a system designed to receive, convert, and dispatch HL7 and DICOM messages, integrating all of your healthcare applications, including imaging modalities, PACS systems and archives. The system supports both DICOM Modality Worklist and DICOM Modality Performed Procedure Step. Unlike other integration engines, UltraBROKER can be used as a stand-alone solution, or it can be integrated, to create new applications or extend your current technologies.



UltraGATEWAY

The Next Generation of DICOM Image Routing!

UltraGATEWAY is a highly configurable router which distributes images and information across point-to-point connections within the facility as well as any remote location for diagnostic or consultation purposes. The flexible software is used widely in the market today with key features every organization needs for clean, efficient and effective workflow.



UltraBURN

The Next Generation of DICOM Media Management!

UltraBURN allows users to burn and print DICOM media from anywhere on the network, or simple rules can automate the process. Easy to use web interface allows the user to customize, configure and control the UltraBURN system.



UltraARCHIVE+

UltraARCHIVE+ provides the ultimate archiving software solution designed to manage the storage of DICOM information for any hospital and imaging center regardless of size. UltraARCHIVE+ delivers exceptional speed, quality, security, simplicity, flexibility, redundancy, and scalability. Start small, and easily expand as your practice grows!



UltraSCAN

UltraSCAN is the perfect choice to help optimize your workflow! Scan critical information, such as Physicians Orders, Lab Reports, EKG's, and other pertinent paper documentation and attach them to the appropriate DICOM patient studies in just a few simple steps.



Ultra3D

Ultra3D solution provides daily tools which allow the user to speed up diagnosis and facilitate the communication with referring physicians and other specialists. The system can be stand alone or fully integrated within your existing workflow. Ultra3D is based on a unique GPU technology software platform which uses standard PCs and off the shelf video cards.



UltraORTHO

State-of-the-art Digital Tools for Surgical Planning!

Orthopedists can view images, take measurements, and plan surgery, without the use of film. With this web based solution, users can make advanced orthopedic measurements that minimize calculations and human error while facilitating quick and accurate surgical planning.




Image Management



UltraWEB

A Complete Image Viewing Solution!

UltraWEB enables your medical professionals to easily share and view images over the Internet and/or Intranet. The Web-based system includes easy image and report viewing and manipulation, plus security and privacy features to facilitate hospital HIPAA compliance.



UltraVIEW

Scalable Diagnostic Viewing Solutions!

UltraVIEW is a complete and cost-effective stand-alone diagnostic viewing system. The system includes: An Intuitive GUI, Comprehensive Workflow, Orthopedic tools and overlay templates, a Full suite of Advanced Visualization Tools and Customizable Hanging Protocols.



UltraREPORTER

A Complete Digital Imaging Workflow and Reporting solution!

The UltraRAD PACS WORKFLOW and REPORTING solution provides customers with a complete image management and reporting solution based on StructuRad's StudyManager and ReportNow. This integrated solution provides a consolidated worklist that orchestrates the workflow for both image viewing and reporting.



UltraAUDIT

UltraAUDIT HIPAA Compliance Manager

The UltraAUDIT solution enables any imaging facility, Hospital, Clinic, Teleradiology service, or other DICOM-based sites to monitor and track study distribution and audit study access of Protected Healthcare Information. The solution assists in HIPAA compliance and supplies tools for a facility to track studies and their usage.