




One Company. One Product. Multiple Solutions.


Upgrade Your Radiographic System to the Finest DR Solution on the Market Today.



Upgrade From Your Film or CR Based Radiographic System Today and Enjoy These benefits:

- Optimum Image Quality from the 9.3 Megapixel, 16 bit A/D Full-Field Imaging Detector.
- Full integration with CPI radiographic generators*
- Cassette Sized Portable Detector for easy and cost effective retrofit. Works with your existing equipment without expensive modifications.
- APR Tuned Image Processing and Generator Techniques.
- Highly automated single console integration that can begin and complete an examination in as few as 4 clicks.


Two Panel Retrofit Solution



Fixed Varian 17 by 17 Detector for Upright Exams.

Portable Varian 14 by 17 Detector for Tabletop and Table Bucky Exams

Single Panel Retrofit Solution



Portable Varian 14 by 17 Detector All Upright, Tabletop and Table Bucky Exams

Note: CPI neither provides nor endorses any positioner, accessory or peripheral other than radiographic and R&F generators and Radiographic digital imaging systems. Images used above are for illustration of commonly used radiographic solutions.

* Contact CPI for details.



45 River Drive
Georgetown, ON
Canada
L7G 2J4
Tel. (905) 877-0161
Fax. (905) 877-5327



www.cpii.com